



# Evaluation of Pyrrolidonyl Arylamidase Activity in *Staphylococcus delphini*

Samantha T. Compton,<sup>a</sup> Stephen A. Kania,<sup>a</sup> Amy E. Robertson,<sup>b</sup>

Sara D. Lawhon,<sup>c</sup> Stephen G. Jenkins,<sup>d</sup> Lars F. Westblade,<sup>d</sup> David A. Bemis<sup>a</sup>

Department of Biomedical and Diagnostic Sciences, University of Tennessee College of Veterinary Medicine, Knoxville, Tennessee, USAa; NewYork–Presbyterian/Weill Cornell Medical Center, New York, New York, USAb; Department of Veterinary Pathobiology, Texas A&M University College of Veterinary Medicine, College Station, Texas, USAc; Department of Pathology and Laboratory Medicine, Weill Cornell Medical College, New York, New York, USAd

**ABSTRACT** Clinical reference textbooks lack data for pyrrolidonyl arylamidase (PYR) activity in *Staphylococcus delphini*. This study evaluated PYR activities of 21 *S. delphini* strains by reference broth, rapid disc, and rapid slide methods. Species and subgroup identifications were confirmed by nucleic acid-based methods and included nine group A and 12 group B strains. Testing by rapid PYR methods with products from four manufacturers was performed at two testing locations, and, with the exception of one strain tested at one location using reagents from one manufacturer, each *S. delphini* strain tested positive for PYR activity. Therefore, PYR may be a useful single-test adjunct for distinguishing *Staphylococcus aureus* from *S. delphini* and other members of the *Staphylococcus intermedius* group.

**KEYWORDS** Staphylococcus delphini, Staphylococcus intermedius group, pyrrolidonyl arylamidase

yrrolidonyl arylamidase (PYR), also known as pyrrolidonyl aminopeptidase, is a bacterial enzyme that hydrolyzes L-pyroglutamic acid-β-naphthylamide to produce β-naphthylamine, which combines with N,N-dimethylaminocinnamaldehyde reagent to form a red color (1). Rapid disc or slide tests for PYR activity are used to differentiate Enterococcus species from Streptococcus species (1) but can also be used to presumptively differentiate Staphylococcus aureus from phenotypically similar staphylococci (2, 3). One such group of staphylococci is the Staphylococcus intermedius group (SIG), which includes the species Staphylococcus delphini, Staphylococcus intermedius, and Staphylococcus pseudintermedius (4). S. delphini has been further subdivided into two distinct phylogenetic clades referred to as groups A and B (5). PYR test results were not provided in the original description of S. delphini (6); however, a subsequent report indicated that 100% of 22 S. delphini isolates (17 group A and 5 group B) were PYR positive according to a commercial multibiochemical test strip (Rapid ID 32 Staph [bioMérieux, Durham, NC, USA]) (5). As highlighted in Table 1, the PYR activity of S. delphini is listed in some clinical reference texts as "not determined" or "not available," and in other references either S. delphini or the PYR test result is excluded from Staphylococcus species identification tables (7-10). In comparison, the PYR activities of S. intermedius and S. pseudintermedius are clearly documented as positive (7-9). Furthermore, rather than using multibiochemical test strips, PYR activity is typically evaluated in the routine clinical microbiology setting with rapid PYR disc or slide kits, which can be used in less than 5 min. Therefore, the purpose of this study was to determine the PYR activities of *S. delphini* strains using the reference standard broth method (3) to complete an important piece of "missing" biochemical data and to evaluate the performance of four rapid PYR test kits produced by different manufacturers to assist

Received 10 October 2016 Returned for modification 1 November 2016 Accepted 19 December 2016

**Accepted manuscript posted online** 21 December 2016

Citation Compton ST, Kania SA, Robertson AE, Lawhon SD, Jenkins SG, Westblade LF, Bemis DA. 2017. Evaluation of pyrrolidonyl arylamidase activity in *Staphylococcus delphini*. J Clin Microbiol 55:859–864. https://doi.org/ 10.1128/JCM.02076-16.

**Editor** Karen C. Carroll, The Johns Hopkins University School of Medicine

**Copyright** © 2017 American Society for Microbiology. All Rights Reserved.

Address correspondence to Lars F. Westblade, law9067@med.cornell.edu, or David A. Bemis, bemis@utk.edu.

Compton et al. Journal of Clinical Microbiology

TABLE 1 Documented PYR activities of S. aureus and members of the S. intermedius group represented in clinical reference texts


	Documented PYR activity <sup>a</sup>				
Clinical reference text	S. aureus	S. delphini	S. intermedius	S. pseudintermedius	
Manual of Clinical Microbiology, 11th ed (7)	NEG	Not determined	POS	POS	
Koneman's Color Atlas and Textbook of Diagnostic Microbiology, 7th ed (8)	NEG	Not available	POS	POS	
Bailey and Scott's Diagnostic Microbiology, 13th ed (9)	NEG	Excluded	POS	POS	
Clinical Veterinary Microbiology, 2nd ed (10)	Excluded from identification table				

<sup>&</sup>lt;sup>a</sup>PYR, pyrrolidonyl arylamidase; NEG, negative; POS, positive.

clinical microbiologists in the rapid differentiation of *S. aureus* from the phenotypically similar coagulase-positive SIG members.

### RESULTS

Each PYR test kit reacted as expected with the control organisms (*Streptococcus agalactiae*, negative; *Enterococcus faecalis*, positive). Only one manufacturer (Key Scientific Products) included expected reactions for *Staphylococcus* species in its instructions. Nonetheless, the selected *S. aureus* (negative) and *S. pseudintermedius* (positive) controls performed consistently with complete agreement for each test kit, at each test point, and, in the case of the rapid PYR disc/slide tests, between laboratories. Positive PYR reactions observed in the reference broth test were indicated with a brilliant red-fuchsia color, while negative reactions were indicated by yellow or weakly orange color (Fig. 1). Positive PYR reactions observed with the rapid disc and rapid slide tests were of a similar red-pink color but were generally weaker than those observed with the reference broth, and in these tests, *Staphylococcus* species (*S. pseudintermedius* and *S. delphini*) produced weaker reactions than the control *Enterococcus* species. Results of PYR tests for the 21 *S. delphini* and control strains are summarized in Table 2. Our results show excellent agreement among the test systems (broth, rapid disc, and rapid slide),



**FIG 1** Observed broth PYR reactions for the control and *S. delphini* DSM 20771<sup>T</sup> strains. 1, *Enterococcus faecalis* ATCC 29212; 2, *Streptococcus agalactiae* ATCC 12386; 3, *Staphylococcus aureus* ATCC 25923; 4, *Staphylococcus pseudintermedius* LMG 22219<sup>T</sup>; 5, *Staphylococcus delphini* DSM 20771<sup>T</sup>.

TABLE 2 Results of PYR broth and rapid PYR disc/slide testing at the two testing sites

	No. of strains tested	Reference PYR broth <sup>b</sup> (% of positive tests)	Rapid PYR disc/slide test results by manufacturer <sup>c</sup> (% of positive tests)			
Organism <sup>a</sup>			A	В	С	D
E. faecalis ATCC 29212	1	100	100	100	100	100
S. agalactiae ATCC 12386	1	0	0	0	0	0
S. aureus ATCC 25923	1	0	0	0	0	0
S. pseudintermedius LMG 22219 <sup>™</sup>	1	100	100	100	100	100
S. delphini (group A)	9	100	100	100	100	94.4 <sup>d</sup>
S. delphini (group B)	12	100	100	100	100	100
S. delphini (groups A and B)	21	100	100	100	100	97.6e

aLMG, Collection of the Laboratorium voor Microbiologie en Microbielle Genetica.

manufacturers, and laboratories, with all *S. delphini* strains testing positive for PYR activity except for one strain that tested negative with one rapid disc test kit in one laboratory. Interlaboratory agreement, based on total number of rapid disc and rapid slide tests (42 rapid disc/slide test results in total for both testing sites) performed on all *S. delphini* strains was 97.6% (41/42 rapid disc/slide test results) (Table 2).

#### **DISCUSSION**

Of the tested S. delphini isolates in this study, 100% tested positive using the reference PYR broth method, while >95% of the rapid disc/slide test results for the 21 S. delphini isolates were positive. The ability to accurately differentiate Staphylococcus aureus from other phenotypically similar coaqulase-positive staphylococci is important for accurate interpretation of antimicrobial susceptibility test results that may influence therapy selection. This has been demonstrated most clearly with respect to the detection of methicillin susceptibility in S. aureus versus S. pseudintermedius (11). While some reference and larger hospital-based microbiology laboratories have turned to DNA sequencingor mass spectrometry-based technologies to rapidly and accurately identify Staphylococcus species (12-14), many laboratories will continue to utilize single rapid tests and short test sets for the identification of S. aureus and SIG members. The use of PYR testing to distinguish some coagulase-positive Staphylococcus species from S. aureus is well established (7, 8, 9; see also http://www.aavmc.org/data/files/case-study/burnham%20-%20staph %20pseudintermedius%20-%20student%20materials.pdf). The rapid PYR disc methods took <5 min to complete, while the reference broth method required 4 h. It should be noted that strict adherence to manufacturer recommendations is required for all PYR-based tests. It is especially important that PYR test reactions be read promptly after adding the development reagent (N,N-dimethylaminocinnamaldehyde), since some negative reactions may appear positive after increased time. Yellow-, salmon-, or orange-colored reactions are generally considered by all manufacturers to be negative. It is recommended that isolates exhibiting pale-pink or weak reactions be retested either with a longer substrate incubation of 5 min, directly with reagent alone to determine if the reactions are specific to the development reagent rather than the substrate, or using the reference tube method.

The purpose of this study was to complete "missing" biochemical data for *S. delphini* to facilitate the rapid differentiation of *S. aureus* and members of the SIG isolated from animal and human clinical specimens. While *S. pseudintermedius* is increasingly recognized as a potential human pathogen (11, 15), *S. delphini* and *S. intermedius* are primarily animal pathogens. To the best of the authors' knowledge, neither *S. delphini* nor *S. intermedius sensu stricto* has been reported as a cause of human infection; likewise, methicillin resistance has not been reported in these species. Each of the *S. delphini* isolates used in this study was presumed to be susceptible to oxacillin by the

<sup>&</sup>lt;sup>b</sup>Thermo Fisher Scientific (Remel); tested only at the University of Tennessee.

cManufacturers are A, Thermo Fisher Scientific (Remel); B, Becton, Dickinson and Company; C, Hardy Diagnostics; and D, Key Scientific Products. All rapid test kits were tested once in each of two laboratories.

<sup>&</sup>lt;sup>d</sup>One *S. delphini* group A isolate, *S. delphini* MI 09-2894, tested negative using the Key Scientific Products PYR disc kit at Weill Cornell Medical College. Of a total of 18 tests results available from both testing sites for the *S. delphini* group A isolates (9 results from each site), 94.4% (17/18) were positive.

eOf 42 test results from both testing sites available for the S. delphini group A and B isolates combined (21 results from each testing site), 97.6% (41/42) were positive.

Compton et al. Journal of Clinical Microbiology

TABLE 3 Staphylococcus delphini strains used in this study

Strain	Group	Host species	Geographic origin	Reference(s) or source
DSM 20771 <sup>Ta</sup>	Α	Dolphin	Italy	6
VPP 12 D	Α	Horse	Switzerland	This study
MI 09-8445	Α	Sea otter	USA	This study
MI 09-2894	Α	Raccoon	USA	This study
MI 16-1129	Α	Ferret	USA	This study
19-039	Α	Avian	USA	11
52-006	Α	Horse	USA	This study
56-021	Α	Llama	USA	This study
62-031	Α	Horse	USA	This study
8086	В	Horse	UK	20
AV8047	В	Pigeon	Japan	20, 21
H4A	В	Horse	Japan	5, 20
HT2003-0674	В	Camel	France	20, 22
P27B	В	Pigeon	Japan	5, 20
P50	В	Pigeon	Japan	5, 20
26-037	В	Horse	USA	11
33-090	В	Horse	USA	11
35-016	В	Horse	USA	11
49-046	В	Horse	USA	This study
53-033	В	Horse	USA	This study
56-039	В	Horse	USA	This study

<sup>&</sup>lt;sup>a</sup>DSM 20771<sup>T</sup>, Deutsch Sammlung von Mikroorganismen und Zellkulturen GmbH, S. delphini type strain.

standard disk diffusion method (range of zone of inhibition diameters, 22 to 28 mm; arithmetic mean, 26 mm) according to interpretive criteria used for *S. pseudintermedius* (11), and the *mecA* gene was not detected by PCR (data not shown). Perhaps due to misidentifications as *S. aureus* (12, 16) and the lack of routine species differentiation within the SIG in some laboratories, the true epidemiology and distribution of *S. delphini* infections in the human population may not be known. *S. delphini* has been isolated from many different animal species (17–19), including some, such as horses, that may have close human contact. Therefore, methods that can be used for quick presumptive recognition of all members of the SIG (including *S. delphini*), such as PYR, will be very useful in medical and veterinary clinical microbiology laboratories.

# **MATERIALS AND METHODS**

**Bacterial strains.** We used four control strains: Enterococcus faecalis strain ATCC 29212, Streptococcus agalactiae strain ATCC 12386, Staphylococcus aureus strain ATCC 25923 (American Type Culture Collection, Manassas, VA, USA), and Staphylococcus pseudintermedius strain LMG 22219 $^{\rm T}$  (received from Freddy Haesebrouck, Ghent University, Ghent, Belgium). The Staphylococcus delphini strains used in this study are summarized in Table 3. The S. delphini strain collection included seven previously published strains, including the S. delphini type strain (S. delphini DSM 20771 $^{\rm T}$ ) (received from Vincent Perreten, University of Bern, Bern, Switzerland), and 14 clinical strains identified in this study. All isolates were cultured on tryptic soy agar with 5% sheep blood (TSAB) (Becton, Dickinson and Company, Franklin Lakes, NJ, USA). Initially, isolates (21 test strains and four control strains) were subcultured from long-term stocks and incubated at 35  $\pm$  2°C in 5% carbon dioxide (CO<sub>2</sub>) for 18 to 24 h. Subsequently, these cultures were subcultured and incubated at 35  $\pm$  2°C in 5% CO<sub>2</sub> for 18 to 24 h prior to PYR testing.

**Molecular identification of Staphylococcus delphini isolates.** Bacterial cell lysates were obtained from colonies grown on TSAB following incubation at  $35 \pm 2^{\circ}\text{C}$  in 5% CO $_2$  for 18 to 24 h. A single isolated colony was suspended in 0.5 ml of a 1:1 mix (vol/vol) containing Tris-EDTA buffer (10 mM Tris, 1 mM EDTA [pH 8.0]) and 0.1-mm zirconium beads (Biospec Products, Bartlesville, OK, USA) and subsequently vortexed to disrupt bacterial cells. The species and subgroup identification of *S. delphini* strains were confirmed by thermonuclease (nuc) gene PCR using previously described primers (23). The total reaction mixture for each PCR was  $25 \,\mu$ l, containing  $2.5 \,\mu$ l of supernatant from the bacterial cell lysate,  $1 \,\mu$ l of each primer (final concentration, 10 pmol),  $8 \,\mu$ l of water, and  $12.5 \,\mu$ l of  $2\times$  PCR solution with a final concentration of  $0.625 \,\text{U}$  rTaq polymerase,  $0.2 \,\text{mM}$  deoxynucleoside triphosphate mixture, and  $1\times$  reaction buffer (Premix Taq, TaKaRa Bio Inc., Kyoto, Japan). The thermocycler parameters were 1 cycle at 95°C for 1 min 30 s, 30 cycles at  $56^{\circ}$ C for 30 s,  $72^{\circ}$ C for 2 min 30 s and  $94^{\circ}$ C for 1 min, 1 cycle at  $50^{\circ}$ C for 2 min, and 1 cycle at  $72^{\circ}$ C for 5 min, followed by holding at  $4^{\circ}$ C. Resultant PCR products were resolved by agarose gel electrophoresis (1.0%) and visualized by staining with ethidium bromide (0.5  $\mu$ g/ml).

**Antimicrobial susceptibility testing.** Each *S. delphini* isolate was tested for oxacillin susceptibility using the standard disk diffusion method as recommended by the Clinical and Laboratory Standards Institute (CLSI) (24). Diameters of the zone of growth inhibition were recorded and, because no interpretive breakpoints that apply specifically to this species exist, comparisons were made to the

current breakpoints used for *S. pseudintermedius* (11, 24). Bacterial cell lysates, obtained as described above, were subjected to PCR for the *mecA* gene using previously described primers, thermocycling conditions, and agarose gel electrophoresis method (25).

PYR testing. Four control strains were used to validate products each time 5. delphini strains were tested. Reference broth testing was performed using Remel PYR broth (catalog number 062085; Thermo Fisher Scientific [Remel], Lenexa, KS, USA) at one location (University of Tennessee). Each isolate was tested only once, and testing was performed and interpreted by a single individual according to manufacturer instructions. PYR broth was inoculated with 4 colonies from an overnight TSAB plate culture and incubated in an aerobic atmosphere for 4 h at 37°C. After incubation, 1 drop of PYR reagent (catalog number R21258) was added to the tube. A red (positive) color reaction was recorded after 1 min and checked for significant change at 2 min. The following rapid disc or rapid slide tests were tested: manufacturer A, PYR disc with reagent (catalog number R30854301; Thermo Fisher Scientific [Remel]); manufacturer B, BBL DrySlide PYR kit (catalog number 231747; Becton, Dickinson and Company); manufacturer C, PYR test kit and reagent (catalog number Z175; Hardy Diagnostics, Santa Maria, CA, USA); and manufacturer D, PYR discs (with PEP) (catalog number K1538B; Key Scientific Products, Stamford, TX, USA). For the rapid disc or rapid slide tests, testing was conducted independently at two locations (University of Tennessee and Weill Cornell Medical College) using kits from identical lots. At each site, isolates were assayed only once by each of the different disc or slide tests. Furthermore, testing was performed and interpreted by a single individual at each site according to manufacturer instructions. The sites were not blinded to the PYR activity of the control strains for the rapid PYR tests, as this served as quality control prior to testing. However, all other testing factors were blinded. Discs and slides were moistened with a 10-µl loopful of distilled water, the incubation time of the organism on discs and slides for all kits prior to addition of reagent was 2 min (including Key Scientific Products, which indicates that the incubation time can be between 2 and 5 min), and the color reaction was recorded 1 min after addition of the reagent for all kits tested (including Key Scientific Products, which instructs at least 1 min but not more than 2 min before recording).

#### **ACKNOWLEDGMENTS**

We thank the representatives of Thermo Fisher Scientific (Remel), Becton, Dickinson and Company, Hardy Diagnostics, and Key Diagnostics for generously providing samples of PYR test products used in this study. We are also grateful to Rupal Brahmbhatt, Mary Jean Bryant, Rebekah Jones, and Brian Johnson for technical assistance, Phil Snow for photographic assistance, and Vincent Perreten and Freddy Haesebrouck for providing reference strains.

S.T.C. was supported by a fellowship from the University of Tennessee College of Veterinary Medicine, Center of Excellence in Livestock Diseases and Human Health, Student Summer Research Program. This research received no specific grant from any funding agency in the public, commercial, or not-for-profit sectors.

## **REFERENCES**

- Facklam RR, Thacker LG, Fox B, Eriquez L. 1982. Presumptive identification of streptococci with a new test system. J Clin Microbiol 15:987–990.
- Oberhofer TR. 1986. Value of the L-pyrrolidonyl-β-napthylamine hydrolysis test for identification of select Gram-positive cocci. Diagn Microbiol Infect Dis 4:43–47. https://doi.org/10.1016/0732-8893(86)90055-6.
- Frank KL, del Pozo JL, Patel R. 2008. From clinical microbiology to infection pathogenesis: how daring to be different works for *Staphylo-coccus lugdunensis*. Clin Microbiol Rev 21:111–133. https://doi.org/ 10.1128/CMR.00036-07.
- Devriese LA, Vancanneyt M, Baele M, Vaneechoutte M, De Graef E, Snauwaert C, Cleenwerck I, Dawyndt P, Swings J, Decostere Haesebrouck AF. 2005. Staphylococcus pseudintermedius sp. nov., a coagulase positive species from animals. Int J Syst Evol Microbiol 55:1569–1573. https:// doi.org/10.1099/ijs.0.63413-0.
- Sasaki T, Kikuchi K, Tanaka Y, Takahashi N, Kamata S, Hiramatsu K. 2007. Reclassification of phenotypically identified *Staphylococcus intermedius* strains. J Clin Microbiol 45:2770–2778. https://doi.org/10.1128/JCM .00360-07.
- Varaldo PE, Kilpper-Balz R, Biavasco F, Satta G, Schleifer KH. 1988. Staphylococcus delphini sp nov. a coagulase-positive species isolated from dolphins. Int J Syst Evol Microbiol 38:436–439.
- 7. Becker K, Skov RL, von Eiff C. 2015. *Staphylococcus, Micrococcus*, and other catalase-positive cocci, p 354–382. *In* Jorgensen JH, Pfaller MA, Carroll KC, Funke G, Landry ML, Richter SS, Warnock DW (ed), Manual of clinical microbiology, 11th ed, vol 1. ASM Press, Washington, DC.
- 8. Procop GW, Church DL, Hall GS, Janda WM, Koneman EW, Schreckenberger PC, Woods GL (ed). 2017. Koneman's color atlas and textbook of

- diagnostic microbiology, p 670–732, 7th ed. Wolters Kluwer, Philadelphia, PA.
- Tille P (ed). 2014. Bailey and Scott's diagnostic microbiology, p 232–246, 13th ed. Mosby. St. Louis. MO.
- Markey B, Leonard F, Archambault M, Cullinane A, Maguire D. 2013. Clinical veterinary microbiology, p 105–119, 2nd ed. Mosby, St. Louis, MO.
- Wu MT, Burnham C, Westblade LF, DienBard J, Lawhon SD, Wallace MA, Stanley T, Burd E, Hindler J, Humphries RM. 2016. Evaluation of oxacillin and cefoxitin disk and MIC breakpoints for prediction of methicillin resistance in human and veterinary isolates of *Staphylococcus interme*dius group. J Clin Microbiol 54:535–542. https://doi.org/10.1128/ JCM.02864-15.
- Viau R, Hujer AM, Hujer KM, Bonomo RA, Jump RLP. 2015. Are Staphylococcus intermedius infections in humans cases of mistaken identity? A case series and literature review. Open Forum Infect Dis 2:ofv110. https://doi.org/10.1093/ofid/ofv110.
- Decristophoris P, Fasola A, Benagli C, Tonolla M, Petrini O. 2011. Identification of Staphylococcus intermedius group by MALDI-TOF MS. Syst Appl Microbiol 34:45–51. https://doi.org/10.1016/j.syapm.2010.11.004.
- Westblade LF, van Belkum A, Grundhoff A, Weinstock GM, Palmer EG, Pallen MJ, Dunne WM, Jr. 2016. Role of clinicogenomics in infectious disease diagnostics and public health microbiology. J Clin Microbiol 54:1686–1693. https://doi.org/10.1128/JCM.02664-15.
- Frank MG, Keniston A, Madinger N, Price C, Bessesen MT. 2015. Staphylococcus intermedius group infections in humans: report of four cases and a literature review. JMM Case Rep 2 https://doi.org/10.1099/jmmcr .0.000066.

Compton et al. Journal of Clinical Microbiology

- Börjesson S, Gomez-Sanz E, Ekström K, Torres C, Grönlund U. 2015. Staphylococcus pseudintermedius can be misdiagnosed as Staphylococcus aureus in humans with dog bite wounds. Eur J Clin Microbiol Infect Dis 34:839 – 844. https://doi.org/10.1007/s10096-014-2300-y.
- Ben Zakour NL, Beatson SA, van den Broek AH, Thoday KL, Fitzgerald JR. 2012. Comparative genomics of the *Staphylococcus intermedius* group of animal pathogens. Front Cell Infect Microbiol 2:44. https://doi.org/ 10.3389/fcimb.2012.00044.
- Guardabassi L, Schmidt KR, Petersen TS, Espinosa-Gongora C, Moodley A, Agersø Y, Olsen JE. 2012. Mustelidae are natural hosts of Staphylococcus delphini group A. Vet Microbiol 159:351–353.
- Stull JW, Slavic D, Rousseau J, Weese JS. 2014. Staphylococcus delphini and methicillin-resistant S. pseudintermedius in horses, Canada. Emerg Infect Dis 20:485–448. https://doi.org/10.3201/eid2003.130139.
- Bannoehr J, Ben Zakour NL, Waller AS, Guardabassi L, Thoday KL, van den Brock AHM, Fitzgerald JR. 2007. Population genetic structure of the Staphylococcus intermedius group: insights into agr diversification and the emergence of methicillin-resistant strains. J Bacteriol 189: 8685–8692. https://doi.org/10.1128/JB.01150-07.
- 21. Futagawa-Saito K, Sugiyama T, Karube S, Sakurai N, Ba-Thein W, Fuku-yasu T. 2004. Prevalence and characterization of leukotoxin-producing

- Staphylococcus intermedius in isolates from dogs and pigeons. J Clin Microbiol 42:5324–5326. https://doi.org/10.1128/JCM.42.11.5324-5326.2004.
- Bes M, Saidi Slim I, Becharnia F, Meugnier H, Vandenesch F, Etienne J, Freney J. 2002. Population diversity of *Staphylococcus intermedius* isolates from various host species: typing by 16S-23S intergenic ribosomal DNA spacer polymorphism analysis. J Clin Microbiol 40:2275–2277. https://doi.org/10.1128/JCM.40.6.2275-2277.2002.
- Sasaki T, Tsubakishita S, Tanaka Y, Sakusabe A, Ohtsuka M, Hirotaki S, Kawakami T, Fukata T, Hiramatsu K. 2010. Multiplex-PCR method for species identification of coagulase-positive staphylococci. J Clin Microbiol 48:765–769. https://doi.org/10.1128/JCM.01232-09.
- 24. Clinical and Laboratory Standards Institute. 2015. Performance standard for antimicrobial disk and dilution susceptibility tests for bacteria isolated from animals (approved standard); 4th edition and supplement. VET01-A4 and VET01S-Ed3. Clinical and Laboratory Standards Institute, Wayne, PA.
- Bemis DA, Jones RD, Frank LA, Kania SA. 2009. Evaluation of susceptibility test breakpoints used to predict mecA-mediated resistance in Staphylococcus pseudintermedius isolated from dogs. J Vet Diagn Invest 21:53–58. https://doi.org/10.1177/104063870902100108.